



Saliva Absorbing Mask for Difficulty Swallowing Patients


PROBLEM:

Patients with ALS, stroke, and Parkinson's disease often experience the **problem of drooling**, which can cause skin discomfort and psychological embarrassment, affecting their **dignity** and willingness to participate in outdoor activities or social events.



METHOD:

Mask with "U" Pocket + Drool Pad




"U" pocket design
for drool pad storage



Easy to replace
drool pads or tissues



Absorbent drool pad
material: cotton



Waterproof to
avoid saliva leakage

BENEFITS:

- Absorb saliva and self-expression
- Easy to replace drool pads or tissues.
- "U" opening pocket: **more contact area** for drool pad to touch the skin.
- The mask and the drool pad are **comfortable & washable** to reuse.
- Improves the patient's **willingness** to participate in social activities.



Shaping the jawline
and various appearance